

TELLUS

XRD series

Multi-function solution with exceptional materials research & analysis speed, reliability and reproducibility X-ray diffractometer that meets the highest requirements for this class of instruments and even exceeds expectations

DESIGNED FOR EDUCATION, RESEARCH AND INDUSTRY

 **LINEV SYSTEMS**
United for Innovation




From Classroom to Advanced
Scientific Research:
Precision You Can Trust



Tailored Solutions for Every Stage of Innovation


The TELLUS XRD Series by LINEV Systems offers cutting-edge X-ray Diffraction solutions, combining precision, efficiency and versatility in a compact design

Developed for education, advanced research and industrial production the TELLUS series deliver unparalleled results across diverse applications



TELLUS

XRD series



TELLUS BASIC

PERFECT FOR EDUCATION

Simplifies XRD concepts for students with an affordable and easy-to-use design

ROUTINE POWDER ANALYSIS

Ideal for entry-level industrial labs

BENCHTOP

Best for cost-effective, routine QC but lacks advanced features

TELLUS PRO

PRECISION FOR RESEARCH AND INDUSTRY

Features higher accuracy and extended functionality beyond the powder analysis

VERSATILITY

Supports applications like thin-film analysis, residual stress and texture measurements

MULTIPURPOSE

Balances flexibility and precision for R&D labs

TELLUS HR

HIGH RESOLUTION EXPERIMENTS

Advanced high-resolution analysis with cutting-edge technology

RESEARCH GRADE

High-precision measurements of wide range of materials and analysis: epitaxial structures, texture, residual stress etc.

MULTIFUNCTIONALITY & ACCURACY

Unmatched resolution and customization

WHICH TELLUS IS RIGHT FOR YOU?

TELLUS
XRD series

TELLUS BASIC

Essential benchtop XRD for education and routine powder analysis in quality control labs

TELLUS
XRD series


TELLUS PRO

Versatile research-grade diffractometer with thin-film, texture and stress analysis capabilities

TELLUS
XRD series

TELLUS HR

High-resolution XRD system for advanced studies of epitaxial films, residual stress and complex materials



GET IN TOUCH



linevsystems.com
info@linevsystems.com

Start Your Journey with the

TELLUS
XRD series